



**RELAX**  
Replacement Bag

Thank you for your purchase of Rainshower's Replacement Bag for the Relax Bath Dechlorinator. The Relax Replacement Bag is designed to treat the chlorine in your bath water for 200 baths or one year, whichever occurs first.

**NOTE:** The Relax Replacement Bag may also be used with the Bath-2000 & Bath-3000 bath dechlorinator.s

## INSTRUCTIONS

### Remove the Old Media Bag

1. Locate the four locking tabs located around the middle of the Relax.
2. Press your thumb above one of the locking tabs and press lightly. While the housing is pressed in slightly, gently rotate the top and bottom halves of the Relax housing in opposite directions and separate the halves.
3. Pull the hanging cord down through the hole in one housing half. Remove the cord and old media bag from the housing and dispose in the trash.

### Install the New Media Bag

1. Thread the hanging cord of the new replacement bag up through the hole in one half of the housing and pull it through, allowing the media bag to fit within both halves of the Relax housing.
2. Align one of the locking tabs from one half of the housing with the inside groove on the opposite housing half. This should align the other three locking tabs with their matching grooves.
3. Gently rotate the housing halves so that all four locking tabs lock into place with their matching locking slots.

7251 Cathedral Rock Dr., Las Vegas, NV 89128  
702.483.4410 • sales@rainshower.com

[www.rainshower.com](http://www.rainshower.com)

Patent # US 6,245,230 • F545-RS / 0120